


AD-P-9020

AUTODROP GANTRY ECO

The **perfect entry** into the world of microdispensing



• ADVANTAGES

- dispensing from pico- to microliter
- software controlled xyz-axis system
- wide viscosity range
- easy refilling and cleaning
- no disposable parts -> no follow on costs
- high repetition and positioning accuracy

• EXTEND IT WITH




- dispenser heads or pipettes
- holder for well plates
- holder for glass slides
- vacuum plate, also heated
- humidifier
- customized substrate plates
- HEPA filter
- your ideas ... contact us.

• PERFECT FOR APPLICATIONS LIKE

- R&D purpose
- small batch productions

POWER PROGRESS ACROSS A RANGE OF APPLICATIONS

- drives up to 2x dispenser heads or pipettes
- travel range of 210 mm x 210 mm
- contact-less dispensing in the pico- to microliter range
- perfect and compact size for lab benches
- ready-to-use system including intuitive software
- positioning accuracy $\pm 25 \mu\text{m}$
- repetition accuracy $\pm 10 \mu\text{m}$
- wide range of optional items



AUTODROP GANTRY ECO

TECHNICAL DATA OF THE AD-P-9020

Electronic control unit	BLDC motors, ball screws, x- and z-axis each 50 W, y-axis 2x 50 W
Travel range	x-axis 210 mm, y-axis 210 mm, z-axis 110 mm
Positioning accuracy	± 25 µm each axis
Repetition accuracy	± 10 µm each axis
Resolution	0.1 µm
Acceleration	maximum 500 mm/s ²
Velocity	maximum 75 mm/s
Additional possible load of positioning unit	maximum 2.5 kg
Dimension (including feet)// Weight	W: 580 mm / D: 580 mm / H: 790 mm // ~62 kg

AUTODROP SOFTWARE

All control electronics are operated via the user interface of the included Autodrop Software


Short description Autodrop Software:

- Setting and saving of dispensing parameter
- Macroscript for programming automatic processes
- Print types: Point to point, in flight, vector based graphics
- Import and edit dxf files for printing curves and other structures (path control)
- Observation of dispensing droplets in flight and determination of drop velocity and drop volume
- Image recognition for substrate position and/or fiducial marks

USE THE SYSTEM WITH


• Autopipettes

Inner nozzle diameter	Standard: 50 µm or 70 µm
Droplet diameter	~70 µm up to 85 µm dispensed with standard pulse, can be reduced up to 50% by using pulse modulation mode
Droplet volume	~180 pl up to 320 pl dispensed with standard pulse, can be reduced up to 50% by using pulse modulation mode
Viscosity range	0.4 up to ~20 mPas



• Micro dispenser heads

Inner nozzle diameter	Standard: 50 µm, 70 µm, 100 µm, 120 µm, 140 µm
Droplet diameter	~65 µm up to 120 µm dispensed with standard pulse, can be reduced up to 50% by using pulse modulation mode
Droplet volume	~140 pl up to 900 pl dispensed with standard pulse, can be reduced up to 50% by using pulse modulation mode
Viscosity range	0.4 up to ~10,000 mPas



EXTENSION POSSIBILITIES

- vacuum tables, also heatable
- holder for well plates
- holder options for glass slides
- laminar flow module
- ionizer
- humidifier

We are there for you. Contact us.

microdrop Technologies GmbH
Tycho-Brahe-Kehre 1
22844 Norderstedt
Germany

Tel. +49 40 | 535383-0
Fax +49 40 | 535383-24
info@microdrop.de
www.microdrop.de
www.microdrop-lifescience.de